

RECYCLE GUIDE



THINK OUTSIDE
THE TRASH

labcon[®]
The Earth Friendly[®] Company

Bring GREEN to your lab bench!


Choosing Labcon products dramatically improves the environmental footprint of your laboratory. We produce the highest quality labware with the lowest environmental impact on the market! Our products are made with 93% clean energy, incorporate sustainable materials, minimize plastic waste, and are designed for maximum recyclability. But don't just take our word for it. Hold us accountable!

Over 100 Labcon products are certified with the ACT Environmental Impact Factor Label to provide a transparent look at our manufacturing process and environmental initiatives. Find out more about our ACT-labeled products at: act.mygreenlab.org/consumables.html



1
PET

Tip Refill Transfer Device



2
HDPE

Tip Refill Base (Recyclable Bioplastic)



4
LDPE

Bags from Tubes & Tips



5
PP

No need to remove label

Tip Refill Platforms + Tube and Tip Racks*

*Rack labels are also #5 PP – leave in place for recycling



7
OTHER

Film* on IntegraPak & Individually Wrapped

*Multi-layer film not currently recyclable



info@labcon.com
www.labcon.com

